



GET IN ON THE DEMAND AND OFFER YOUR PATIENTS THE POC THEY WANT!

Restoring your patients' freedom and control of their day is a big part of what makes the OxyGo® so special. Patients can power an OxyGo by plugging it into a wall outlet, DC outlet in their car, or for hours at a time with each rechargeable battery.

The OxyGo makes jumping in the car to run errands or taking a weekend trips easy. Users can come and go as they please without having to switch equipment. Patients can even take long trips that they never thought possible! The OxyGo is made to follow your patients anywhere... whether it's around town, on the road, or on a cruise.

Weighing in at less than 5 pounds, and with an exceptional battery life, the OxyGo is completely wearable.

- Up to 9.0 hours of battery life with optional double battery
- Only 4.8 lbs. — including the battery!
- No carts or tanks required
- FAA-approved for domestic and international air travel

Give us a call to learn more about how OxyGo can keep your patients going, and keep your costs down!

888-327-7301 | orders@applied-inc.com

See How OxyGo® Outperforms The Competition...



#1400-1000
OxyGo with single cell battery

#1400-1000-16
OxyGo with double cell battery

OxyGo® single cell battery	OxyGo® double cell battery	Respironics SimplyGo Mini® std. battery	Respironics SimplyGo Mini® ext. battery	Inova Labs LifeChoice Activox 4L®
--------------------------------------	--------------------------------------	---------------------------------------------------	---------------------------------------------------	------------------------------------------

	1,050 ml/min	1,050 ml/min	1,000 ml/min	1,000 ml/min	450 ml/min
O₂ Capacity	1,050 ml/min	1,050 ml/min	1,000 ml/min	1,000 ml/min	450 ml/min
FAA Approved	✓	✓	✗	✗	✓
Single Solution	✓	✓	✗	✗	✗
User Replaceable Batteries	✓	✓	✓	✓	✗
Intelligent Delivery Technology: Flow Settings: 1,2,3,4,5	✓	✓	✗	✗	✗
Weight	4.8 lbs. 2.2 kg.	5.8 lbs. 2.6 kg.	5.0 lbs. 2.3 kg.	6.0 lbs. 2.7 kg.	4.8 lbs. 2.0 kg.
Battery Duration	UP TO 4.5 hours	UP TO 9.0 hours	UP TO 4.5 hours	UP TO 9.0 hours	UP TO 8.25 hours
Sound Level at setting 2	39 dB.	39 dB.	43 dB.	43 dB.	44 dB.

View the clinical studies: www.bitly.com/OxyGoClinical